ENGLISH

Step 2: Procedure

Attention: New Procedure utilizes a 200 microliter Pipetter not included in the kit, but sold separately. Catalog no. 0019

1. Remove the cassette from the foil pouch. Label each cassette with the corresponding sample identification.
2. Gently mix the sample. Using a clean pipette tip from the kit and a 200 microliter Pipetter, completely depress the plunger and submerge the pipette tip into the sample.
Gently release the plunger, which will allow the sample to be drawn into the pipette tip. Quickly releasing the plunger may result in air being drawn into the pipette tip which can affect the accuracy of the test result.
3. Remove the dilution vial cap and transfer the sample into the dilution vial. Rinse the pipette by repeatedly drawing up and expressing its contents into the dilution vial.
Recap the dilution vial and invert several times to thoroughly mix the sample with the dilution fluid. **Avoid shaking** the mixture.
4. Using the Pipetter and pipette tip from the previous step, completely depress the plunger and submerge the pipette tip into the dilution mixture. Gently release the plunger to allow the diluted sample to be drawn into the pipette tip.
Completely fill the cassette sample well (marked "S") with the diluted sample; make sure that the cassette is on a level surface.
Allow the test to proceed for 20 minutes before reading the results. **Do not read results after 40 minutes.**

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


Step 3:

Interpretation of Test Results

In a properly functioning test, a line will appear at the Control (C) position of the results window.

Adequate Level is indicated by ONE line at the "C" position.

Inadequate Level is indicated by TWO lines: One at the "T" position and One at the "C" position. Irregardless of the line intensity, a faint line at the "T" position should be interpreted as an **inadequate level** of IgG.

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Note: This improved procedure utilizes a 200 microliter Pipetter which is sold separately (catalog no. 0019). Each kit contains a cassette, disposable pipette tips, dilution vial and instructions.

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For more information please contact us:

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Plasma Calf IgG

Midland Quick Test Kit®

A qualitative test of immunity levels to determine the health and value of newborn calves.
US Patent Nos. 6,245,577 and 6,660,534

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14000 Procedure rev. 0
Copyright © 2004 Midland BioProducts Corporation®
All rights reserved.

The complex block features a large illustration of a cow and calf. A purple starburst in the top left corner says 'New Procedure'. A purple banner on the left side contains the text 'Plasma Calf IgG'. The main title 'Midland Quick Test Kit®' is prominently displayed in the center. Below the title, there are four paragraphs of text in different languages, each describing the test's purpose and providing patent information. At the bottom left, there is vertical text providing contact information and copyright details.